



Power Wheelchair Product Guide

A breakdown of the types of power wheelchairs covered by insurance and the makes and models we currently offer





Table of Contents

Power Mobility Device Categories	03
Lightweight Portable Power Wheelchair	05
Standard Power Wheelchair	06
Heavy Duty Power Wheelchair	07
Very Heavy Duty Power Wheelchair	08
Power Wheelchair w/Tilt	10
Heavy Duty Power Wheelchair w/Tilt	11
Power Wheelchair w/Tilt and Recline	13
Heavy Duty Power Wheelchair w/Tilt & Recline	14
Very Heavy Duty Power Wheelchair w/Tilt and Recline	15

What Category Of Power Mobility Device Is Right For You?

Here is a high level overview of the different categories of power wheelchairs to help you understand which may be most suitable for you. For personalized advice on determining the most suitable power mobility device, consulting a medical professional is strongly recommended. Don't hesitate to give our team a call to learn more about qualifying for a power mobility device through your insurance.



Power Wheelchair

- Individuals who experience limited functional mobility due to weak lower and upper extremity strength, pain, poor endurance, and are at risk of falling
- Individuals who require support to maintain stability and body position during mobility to perform daily activities



Power Wheelchair with Tilt Only (Most Versatile):

- Individuals who require support to change body position to redistribute weight and who are at high risk of developing skin breakdown at the lower back and buttocks due to prolonged seating
- Individuals who are diagnosed with medical condition(s) that include any type of brain disorder, muscle disorder, metabolic disorder, increased muscle tone, muscle weakness, spinal stenosis, stroke, quadriplegia, hemiplegia, paralysis, neuropathy, etc.



Power Wheelchair with Tilt and Recline

- Individuals who require support to change body position to redistribute weight and who are at high risk of developing skin breakdown at the lower back and buttocks due to prolonged seating
- Individuals who require frequent dynamic body repositioning to a wide range of positions for activities such as breathing, swallowing, transferring, pain management, circulation, edema, etc.
- Individuals who are diagnosed with any medical condition(s) that include any type of brain disorder, muscle disorder, metabolic disorder, increased muscle tone, muscle weakness, spinal stenosis, stroke, quadriplegia, hemiplegia, paralysis, neuropathy, use of ventilator, etc.
- When applicable, power assist seat height adjustment feature would aid in non-level transfers and completing activities at higher levels (Shelves, Cabinets, Appliances, etc.)



Power Wheelchair

Lightweight Portable Power Wheelchair

Go Chair MED From PRIDE

- 18" x 17" or 20" x 18", 5-point swivel seat with latch release and 60" lap belt
- Compact turning radius
- Large perforated footplate
- One-hand feather-touch disassembly for easy transport
- The heaviest piece is only 36 lbs

Specifications

Weight Capacity	<u>300 lbs.</u>
Maximum Speed	Up to 3.7 mph
Ground Clearance	1.75" at motor
Seat-To-Ground	22"-23"
Overall Width	25.5" with armrest
Heaviest Piece	36 lbs. (rear section)
When Disassembled	
Range Per Charge	Up to 9.41 miles
Available Colors	●
Add-Ons Covered By Insurance	Oxygen holder
Out Of Pocket Only Add Ons	Cup holder, USB charger, Saddle bag, Weather cover,



Easy disassembly for quick transport




GC3 controller with quick disconnect

Power Wheelchair

Jazzy Elite ES from PRIDE

- Depth-adjustable, limited-recline, fold-flat, high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- Easy-access free-wheel lever
- Specialty seating/cushion options available for those who qualify

Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	1.85" at motors
Overall Width	22.75" - 26.75"
Seat-to-Floor	21.5"-23.5"
Range Per Charge	Up to 15 miles
High-back Seat Sizes	16"W, 18"W, 20"W
Available Colors	
Add-Ons Covered By Insurance	Oxygen holder, Manual elevating leg rest, limb support
Out Of Pocket Only Add Ons	Cup holder, USB charger, Saddle bag, Weather cover, Rear basket



GC3 Controller with quick disconnect




Comfortable high back seat with headrest

Heavy Duty Power Wheelchair

Jazzy Elite HD from PRIDE

- Depth-adjustable, limited-recline, fold-flat, high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Large front wheels absorb and climb, easily tackle life's everyday obstacles.
- Easy-access free-wheel lever

Specifications

Weight Capacity	450 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	3"
Overall Width	24.25" - 28.75"
Seat-to-Floor	20"-22"
Range Per Charge	Up to 14.4 miles
High-back Seat Sizes	20"W, 22"W
Available Colors	
Add-Ons Covered By Insurance	Oxygen holder, Manual elevating leg rest, limb support
Out Of Pocket Only Add Ons	Cup holder, USB charger, saddle bag, Weather cover, Rear basket, Cell phone holder



Depth-adjustable high-back seat



Intuitive, easy to operate controller





Very Heavy Duty Power Wheelchair

Atlantis P710 From Merits

- Super high power in-line motors help rides negotiate 10° uphill ramps with ease
- Four seat posts provide a stable & durable structure for the frame
- Extra-wide seat and extra-high back provide great comfort for bariatric users

Specifications

Weight Capacity	600 lbs.
Maximum Speed	Up to 5 mph
Ground Clearance	3"
Overall Width	30"
Seat-to-Floor	23" - 24"
Range Per Charge	Up to 32 miles
High-back Seat Sizes	24" W
Available Colors	 
Add-Ons Covered By Insurance	Oxygen holder, Manual elevating leg rest, Limb support
Out Of Pocket Only Add Ons	Cup holder, USB charger, Side bag, Rear basket, Cell Phone holder, Weather cover



Impressive 32-mile drive range per charge



Durable & Stable four-post seat base




Power Wheelchair w/ Tilt Only

Power Wheelchair With Tilt

XLR 14 From SHOPRIDER

- Power tilt system enabling independent pressure relief
- Added stability thanks to 6 points of contact with ground
- Self Adjusting 14" Mid-Wheel Drive w/Suspension System
- Manual Semi Reclining Back

Specifications

Weight Capacity	300 lbs.
Maximum Speed	4 mph
Ground Clearance	2.75"
Tilt Degrees	50' of tilt
Total Unit Weight	283 lbs Including Batteries and seating
Battery Options	Group 22NF
Range Per Charge	16 miles
Available Colors	
Add-Ons Covered By Insurance	Swing-away joystick, Oxygen holder, Manual ELR, Arm trough, Stump support
Out Of Pocket Only Add Ons	Lighting kit, Transit tie-downs



Power Tilt Seating System




14" Mid-Wheel Drive w/Suspension

Heavy Duty Power Wheelchair With Tilt

Vision Ultra HD From Merits

- Mid-wheel drive design provides outstanding maneuverability
- Full suspension design offer comfortable ride
- Manual recline, 90°, 95°, 100°, 105° and 110°
- Dual in-line motors for enhanced performance

Specifications

Weight Capacity	450 lbs.
Maximum Speed	5 mph
Ground Clearance	2"
Turning Radius	21"
Overall Height	41"
Overall Length	35.5"
Battery Type	22NF
Range Per Charge	20 mi
Available Colors	
Add-Ons Covered By Insurance	Manual ELR, Oxygen holder, Swing-away joystick, Stump support, Arm trough
Out Of Pocket Only Add Ons	USB charger



Power tilt operated in 0°~55°



Full Suspension Design





Power Wheelchair w/ Tilt & Recline

Power Wheelchair With Tilt & Recline


TDX SP2 From Invacare

- Mid-wheel drive provides tight turning radius for optimal mobility
- 50° of tilt for efficient pressure relief
- 168° of recline with Extended Shear Reduction (ESR) for aid in preventing pressure injuries
- Power Center Mount Elevating Foot Platform for smooth leg elevation
- SureStep suspension and stability lock for smoother transition over obstacles



Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 5.8 mph
Ground Clearance	2.5"
Seat to Floor	17.25", 18.25", 19.25"
Seat Dimensions	16"-22" W x 17"-23" D
Base Width	24"-25.5"
Range Per Charge	Up to 20 miles
Battery Type	Group 22NF
Available Colors	
Add-Ons Covered By Insurance	Oxygen Holder, Individual power ELR, Power seat elevation, & more based on diagnosis
Out Of Pocket Only Add Ons	Cup holder, Lighting kit, Touchscreen joystick, Attendant controller, Transfer handles, Lap tray





Heavy Duty Power Wheelchair Rehab With Tilt & Recline

TDX SP2 HD From Invacare

- Mid-wheel drive provides tight turning radius for optimal mobility
- 50° of tilt for efficient pressure relief
- 168° of recline with Extended Shear Reduction (ESR) for aid in preventing pressure injuries
- Power Center Mount Elevating Foot Platform for smooth leg elevation
- SureStep suspension and stability lock for smoother transition over obstacles

Specifications

Weight Capacity	400 lbs.
Maximum Speed	Up to 5.8 mph
Ground Clearance	2.5"
Seat to Floor	17.75", 18.75", 19.25"
Seat Dimensions	20"-26" W x 17"-23" D
Base Width	25.5"
Range Per Charge	Up to 14 Miles
Battery Type	Group 24
Available Colors	
Add-Ons Covered By Insurance	Oxygen Holder, Individual power ELR, Power seat elevation, & more based on diagnosis
Out Of Pocket Only Add Ons	Cup holder, Lighting kit, Touchscreen joystick, Attendant controller, Transfer handles, Lap tray




Very Heavy Duty Power Wheelchair Rehab With Tilt & Recline

AllTrack HD 3 From AmyLior

- Fully independent 6-wheel interactive suspension with vertical shock absorption tuned for heavier users
- 45° of tilt for efficient pressure relief
- 160° of recline for maintaining posture
- Power center mount elevating foot platform for smooth leg elevation
- Extremely compact turning radius of 20.5"

Specifications

Weight Capacity	550 lbs.
Maximum Speed	Up to 5.4 mph
Ground Clearance	3.5"
Seat to Floor	17"-19"
Seat Dimensions	19"-28" W x 17"-23" D
Base Width	27.4"
Range Per Charge	Up to 14 Miles
Battery Type	Group 24
Available Colors	
Add-Ons Covered By Insurance	Oxygen Holder, Individual power ELR, & more based on diagnosis
Out Of Pocket Only Add Ons	Cup holder, Lighting kit, Touchscreen joystick, Attendant controller, Transfer handles, Lap tray





UniversalMed Supply